

# **Stabilized Alkali Metals for Safer, Scalable Birch Reductions**

Paul F. Vogt, Michael J. Costanzo, Kathryn A. Petersen and Michael Lefenfeld SiGNa Chemistry, Inc., 1 Deer Park Drive, Suite C, Monmouth Junction, NJ 08852

#### Abstract

The Birch reduction is the reduction of aromatic rings with alkali metals (Li, Na, K) dissolved in liquid ammonia. Birch reductions are generally avoided for scale-up due to the complexities of the large-scale process and the hazards associated with the use of the metal and liquid ammonia on scale. SiGNa Chemistry, Inc. has discovered and developed a proprietary method for encapsulating alkali metals in nano-scale porous oxides. The resulting sand-like powders are stable at room temperature and take a significant amount of the danger and associated costs out of using reactive metals. More importantly, the SiGNa materials retain the chemical utility associated with the parent metal. This presentation will highlight SiGNa's technology and how these novel materials can be applied to a range of chemical processes, specifically Birch reductions. A Birch reduction performed using SiGNa material removes the environmental and occupational hazards of alkali metals and does not require liquid ammonia or cryogenic temperatures.

### **The Birch Reduction**



- The Birch reduction<sup>1</sup> is an alternative to hydrogenation that yields cyclohexadienes
- · The reduction is carried out in liquid ammonia with dissolved sodium, lithium or potassium
- · By dissolving the metal in liquid ammonia, a metal cation and a solvated electron are formed. This solvated electron is a powerful reducing agent
- · The reaction mixture contains a stoichiometric amount of an alcohol which acts as a proton dono
- The first step of the mechanism is a one-electron transfer to form a radical anion which is then protonated to yield a cyclohexadienyl radical. This resonance-stabilized allyl radical is converted into a cyclohexadienyl anion by an additional one-electron transfer. Subsequently, the cyclohexadienyl anion is also protonated to yield the 1,4cyclohexadiene

· Choice and handling of materials

· Special design of the equipment

Pungent odor detectable at 5 ppm

Immediate danger to health at 500 ppm

Ammonia3:

# **Issues with Birch Reductions**

- · Handling of alkali metals
- · Handling of ammonia gas
- Cryogenic temperatures

#### Alkali metals2:

- · React Immediately with the environment leading to onerous storage and handling requirements
- Increased risk of supply chain interruption and increased shipping expense

#### A need exists for alkali metals that can be easily handled without a loss in reactivity

# SiGNa's Technology

- Alkali metals encapsulted in nano-structured porous oxides<sup>4</sup>
- Alkali metal silica materials (M-SG); capacity up to 40 wt.%
- Stable in dry air / Safer to handle / easier to transport
- Free-flowing powders

Theoretical Model of Stage I M-SG

- Reducing power of the parent metal/alloy retained
- Non-toxic by-products and waste streams Sodium Silicate

Material	Metals Used	Applications
Stage 0 M-SG	Liquid alloys (K <sub>2</sub> Na, Na <sub>2</sub> K, etc.)	Desulfurization, Polymer Initiation, etc.
Stage I M-SG	Liquid alloys and Na, K, Cs, Rb, etc.	Birch red., deprotections, Wurtz, etc.
Stage II M-SG	Liquid alloys and sodium	Hydrogenation, solvent drying, etc.
Sodium Silicide	Sodium metal	Hydrogen source, fuel cells, etc.





### **Birch Reduction Comparison**



#### Classical Reaction Conditions5.

- · Phenanthrene was dissolved in THF and charged to refluxing ammonia at -33 ° C
- · Lithium metal (3.5 equivalents) was added
- · Reaction quenched with ethanol/water
- · Desired product was isolated in 57% yield
- SiGNa's Reaction Conditions: Phenanthrene and Na-SG(I) (3.5 equivalents)
- are charged to a round bottom flask THF added and the mixture is cooled to 0 °C
- tert-Butanol (1.75 equivalents) is added in one portion
- Desired product was isolated in 60% yield

#### NO CRYOGENIC TEMPERATURES! NO PYROPHORIC REAGENTS! NO LIQUID AMMONIA!

Similar yields and purities were observed with:

- 2.0 to 7.5 equivalents of SG(I) with Na, Na<sub>2</sub>K and K<sub>2</sub>Na as parent metal
- t-butanol, t-amyl alcohol, NH<sub>4</sub>Cl as proton source
- THF, Me-THF, MTBE, Toluene, Cyclohexane, Heptane, or DME as solvent
- Ethylenediamine as an additive

# **Birch Reduction Products**



#### **Synthetic Chemistry Applications**

tetrahydronaphthalene)



## References

- 1. (a) Birch, A. J. J. Chem. Soc. 1944, 430; (b) Birch, A. J. J. Chem. Soc. 1945, 809; (c) Birch, A. J. J. Chem Soc. 1946, 593; (d) Rabideau, P. W. et al. Ore. React. 1992, 42, 1
- 2. Joshi, D. K. et al. Organic Process Research & Development 2005, 9, 997.
- 3. Gangopadhyay, R.K. et al. Process Safety Progress 2008, 27, 15.
- 4. Dye, J. L.; Cram, K. D.; Urbin, S. A.; Redko, M. Y.; Jackson, J. E.; Lefenfeld, M. J. Am. Chem. Soc., 2005, 127, 9338
- 5. P. W. Rabideau, et al. Journal of Organic Chemistry 35(1), 1970, 25.



- · Operational, waste treatment, and safety issues