




The Biotechnology Education Company®



Revised
and
Updated

EDVO-Kit #

221

**Transformation
of *E. coli* with
pGAL™ (blue colony)**

**Storage: See Page 3 for
specific storage instructions**

EXPERIMENT OBJECTIVE:

The objective of this experiment module is to develop an understanding of bacterial transformation by plasmid DNA. This experiment introduces an opportunity to observe an acquired phenotypic trait of the transformed bacterial cells. The presence of blue bacterial colonies visually demonstrates the expression of a specific gene for the Lac⁺ phenotype.

No IPTG used in this experiment.

All components are intended for educational research only.
They are not to be used for diagnostic or drug purposes, nor
administered to or consumed by humans or animals.

The Biotechnology Education Company® • 1-800-EDVOTEK • www.edvotek.com

Table of Contents

	Page
Experiment Components	3
Experiment Requirements	4
Background Information	5
Experiment Procedures	
Experiment Overview	8
Laboratory Safety	10
Experimental Procedures: Transformation of <i>E. coli</i>	11
Experiment Results and Analysis	15
Study Questions	18
Instructor's Guide	
Notes to the Instructor	19
Pre-Lab Preparations	21
Avoiding Common Pitfalls	26
Experiment Results and Analysis	27
Study Questions and Answers	30
Material Safety Data Sheets	31



Important READ ME!

Transformation experiments contain antibiotics which are used for the selection of transformed bacteria. Students who have allergies to antibiotics such as penicillin, ampicillin, kanamycin or tetracycline should not participate in this experiment.

Experiment Components

Contents

**Component
Quantities:**

Experiment # 221
is designed for 10
groups (2-4 students
per group).

Storage: Store components A-G in the refrigerator.

- A Transformation LyphoCells™ (DO NOT FREEZE)
- B Supercoiled pGAL™ (blue colony)
- C Control Buffer (no DNA)
- D Ampicillin
- E X-Gal in solvent (pre-measured)**
- F Cell reconstitution medium
- G Solvent for induction of competency

Storage: Store components listed below at Room temperature

- Bottle ReadyPour™ Luria Broth Agar, sterile (also referred to as ReadyPour medium)
- Bottle Luria Broth Medium for Recovery, sterile (also referred to as Luria Recovery Broth)
- Petri plates
- Plastic microtipped transfer pipets
- inoculating loops (sterile)
- Microtest tubes with attached lids

**EDVOTEK's solvent formulation offers a safer alternative to the usual solvent, N, N-dimethylformamide (DMF).

All components are intended for educational research only.
They are not to be used for diagnostic or drug purposes, nor
administered to or consumed by humans or animals.

None of the experiment components are derived from human
sources.

Requirements

- Automatic Micropipet (5-50 µl) and tips
- Two Water baths (37°C and 42°C)
- Thermometer
- Incubation Oven (34°C and 37°C)
- Pipet pumps or bulbs
- Ice
- Marking pens
- Bunsen burner, hot plate or microwave oven
- Hot gloves

* If a second water bath is not available, water can be heated to 42°C in a beaker. The cells will require this temperature for only a few minutes. Alternatively, 42°C water can be put in a small styrofoam container with a cover. The temperature needs to be held at 42°C.

Online Ordering
now available

www.edvotek.com

Visit our web site for information
about EDVOTEK's complete
line of experiments for biotechnology
and biology education.

Technical Service Department



Mon - Fri
9:00 am to 6:00 pm ET

FAX: (301) 340-0582
web: www.edvotek.com
email: edvotek@aol.com

Please have the following information:

- The experiment number and title
- Kit Lot number on box or tube
- The literature version number
(in lower right corner)
- Approximate purchase date



Bacterial Transformation

Background Information

Bacterial transformation is of central importance in molecular biology. It allows for the introduction of genetically engineered or naturally occurring plasmids in bacterial cells. This makes possible the propagation, genetic expression and isolation of DNA plasmids.

The transformation process involves the uptake of exogenous DNA by cells which results in a newly acquired genetic trait that is stable and heritable. Bacterial cells must be in a particular physiological state before they can be transformed. This state is referred to as competency. Competency can occur naturally in certain species of *Haemophilus* and *Bacillus* when the levels of nutrients and oxygen are low. Competent *Haemophilus* expresses a membrane associated transport complex which binds and transfers certain DNA molecules from the medium into the cell where they are incorporated and their genes are expressed. In nature, the source of external DNA is from other cells.

Most of the current transformation experiments involve *E. coli*. This organism does not enter a stage of competency unless artificially induced. Treatment to achieve competency involves the use of chloride salts, such as calcium chloride, and sudden hot and cold temperature changes. The metal ions and temperature changes affect the structure and permeability of the cell wall and membrane so that DNA molecules can be absorbed by the bacteria. The mechanism of DNA transport in the cell still is not fully understood. Competent *E. coli* cells are fragile and must be treated carefully.

The transformation efficiency is defined by the number of transformants obtained per microgram of DNA. For example, 10 nanograms of DNA were used for a transformation and the cells were allowed to recover in a final volume of 1 ml. One tenth of this volume was plated and produced 100 colonies on a selective agar medium. Therefore, 1000 transformants are present per ml. Keeping in mind that each colony grew from one transformed cell, the efficiency would be $1000/0.01\mu\text{g} = 1 \times 10^5$. Transformation efficiencies of 10^5 to 10^6 are more than sufficient for most subcloning experiments. When the cloning of single copy genes from genomic DNA is done, the required efficiencies are 10^7 to 10^8 .

The determination for transformation efficiency in this case is outlined in Figure 1. Transformation efficiencies gen-

$$\frac{\text{Number of transformants}}{\mu\text{g of DNA}} \times \frac{\text{final vol at recovery (ml)}}{\text{vol plated (ml)}} = \text{Number of transformants per } \mu\text{g}$$

Specific example:

$$\frac{100 \text{ transformants}}{0.01 \mu\text{g}} \times \frac{1 \text{ ml}}{0.1 \text{ ml}} = \frac{100,000 (1 \times 10^5) \text{ transformants}}{\text{per } \mu\text{g}}$$

Figure 1:
Bacterial Transformation Efficiency Calculation

Bacterial Transformation

Background Information

erally range from 1×10^4 to 1×10^7 cells per microgram of DNA. There are special procedures which can produce cells having transformation efficiencies approaching 10^{10} . However, transformation is never 100% efficient. Approximately 1 in every 10,000 cells successfully incorporates plasmid DNA in preparations having average competency. However, there is such a large number of cells in a sample (typically 1×10^9) that only a small fraction needs to be transformed to obtain colonies on a plate. The same volume of recovered cells plated on selective (contains antibiotic) and nonselective agar medium will yield vastly different numbers of cells. The nonselective medium will have many more growing cells that form a bacterial lawn.




Figure 2 :
DNA map of pGAL
Not all restriction enzymes are shown.

Many different plasmids serve as useful tools in molecular biology. One example is the pGAL plasmid, present in multiple copies in specified host *E. coli* host cells. It contains 6751 base pairs and has been cleverly modified by genetic engineering. In the cell, it **does not** integrate into the bacterial chromosome, but replicates autonomously. The pGAL plasmid contains the *E. coli* gene which codes for β -galactosidase. In the presence of artificial galactosides such as 5-Bromo-4 Chloro 3-indolyl- β -D-galactoside (X-Gal), pGAL colonies appear blue when X-Gal is cleaved by β -galactosidase and forms a colored product.

This experiment has been designed to utilize EDVOTEK Transformation LyphoCells™. It also contains the proprietary plasmid, pGAL (Blue Colony), which was engineered by EDVOTEK. Plasmid pGAL carries the complete gene for β -galactosidase. Since the host *E. coli* does not contain a β -galactosidase gene, only cells transformed by the pGAL plasmid will produce the functional β -galactosidase enzyme. Cells that express β -galactosidase will cleave X-Gal and the pGAL transformed colonies will be blue.



Duplication of this document, in conjunction with use of accompanying reagents, is permitted for classroom/laboratory use only. This document, or any part, may not be reproduced or distributed for any other purpose without the written consent of EDVOTEK, Inc. Copyright © 1989, 90, 93, 94, 97, 98, 2000, 2005, EDVOTEK, Inc., all rights reserved
EVT 005077K

Bacterial Transformation

In addition to the expression and cleavage of X-Gal by β -galactosidase, transformation by pGAL is also demonstrated by resistance to ampicillin. *E. coli* host cells used in this experiment are **not** naturally resistant to ampicillin. The plasmid pGAL contains the gene which encodes for β -lactamase that inactivates ampicillin. *E. coli* cells transformed by pGAL will express the resistance gene product β -lactamase as an extracellular enzyme excreted from *E. coli* cells. Once outside the cell, the enzyme diffuses into the surrounding medium and inactivates ampicillin.

With time, small "satellite" colonies may appear around a large blue colony. Cells in the small "satellite" or "feeder" colonies are not resistant to ampicillin and have not been transformed with the pGAL plasmid. They are simply growing in a region of agar where β -lactamase has diffused and inactivated the antibiotic ampicillin. The number of satellite colonies increases if the concentration of ampicillin is low or the plates have incubated for longer times.

Experiment Overview

BEFORE YOU START THE EXPERIMENT

1. Read all instructions before starting the experiment.
2. Write a hypothesis that reflects the experiment and predict experimental outcomes.

EXPERIMENT OBJECTIVE:

The objective of this experiment module is to develop an understanding of the biologic process of bacterial transformation by plasmid DNA. This experiment demonstrates the acquired Lac⁺ phenotypic trait of the transformed bacterial cells as shown by the presence of blue bacterial colonies.

BRIEF DESCRIPTION OF EXPERIMENT:


In this experiment, students will transform host bacterial cells with a plasmid DNA. The transformants acquire antibiotic resistance and exhibit a blue color due to the incorporation and expression of β -galactosidase and ampicillin resistance genes. IPTG is not required since pGAL™ contains the intact β -galactosidase gene. The number of transformants will be counted and the transformation efficiency will be determined.



Duplication of this document, in conjunction with use of accompanying reagents, is permitted for classroom/laboratory use only. This document, or any part, may not be reproduced or distributed for any other purpose without the written consent of EDVOTEK, Inc. Copyright © 1989, 90, 93, 94, 97, 98, 2000, 2005, EDVOTEK, Inc., all rights reserved
EVT 005077K


Experiment Overview

TRANSFORMATION EXPERIMENT FLOW CHART



For optimal results, store covered plates in the upright position after streaking to allow the cell suspension to be absorbed by the agar. After approximately 20 minutes, invert the plates for overnight incubation at 37°C.


Laboratory Safety




Important READ ME!

Transformation experiments contain antibiotics which are used for the selection of transformed bacteria. Students who have allergies to antibiotics such as penicillin, ampicillin, kanamycin or tetracycline should not participate in this experiment.

1. Gloves and goggles should be worn routinely as good laboratory practice.
2. Exercise extreme caution when working with equipment which is used in conjunction with the heating and/or melting of reagents.
3. DO NOT MOUTH PIPET REAGENTS - USE PIPET PUMPS OR BULBS.
4. The *E. coli* bacteria used in this experiment is not considered pathogenic. Although it is rarely associated with any illness in healthy individuals, it is good practice to follow simple safety guidelines in handling and disposal of materials contaminated with bacteria.
5. Properly dispose materials after completing the experiment:
 - A. Wipe down the lab bench with a 10% bleach solution or a laboratory disinfectant.
 - B. All materials, including petri plates, pipets, transfer pipets, loops and tubes, that come in contact with bacteria should be disinfected before disposal in the garbage. Disinfect materials as soon as possible after use in one of the following ways:
 - Autoclave at 121° C for 20 minutes.
Tape several petri plates together and close tube caps before disposal. Collect all contaminated materials in an autoclavable, disposable bag. Seal the bag and place it in a metal tray to prevent any possibility of liquid medium or agar from spilling into the sterilizer chamber.
 - Soak in 10% bleach solution.
Immerse petri plates, open tubes and other contaminated materials into a tub containing a 10% bleach solution. Soak the materials overnight and then discard. Wear gloves and goggles when working with bleach.
6. Wear gloves, and at the end of the experiment, wash hands thoroughly with soap and water.




Duplication of this document, in conjunction with use of accompanying reagents, is permitted for classroom/laboratory use only. This document, or any part, may not be reproduced or distributed for any other purpose without the written consent of EDVOTEK, Inc. Copyright © 1989, 90, 93, 94, 97, 98, 2000, 2005, EDVOTEK, Inc., all rights reserved
EVT 005077K

Transformation of *E. coli*

Guide for Sterile
Calibrated
Transfer Pipet

SETTING UP THE TRANSFORMATION AND CONTROL EXPERIMENT

1. Put your initials or group number on the tubes labeled "pGAL DNA" (contains 25 μ l of plasmid DNA) and "Control Buffer" (contains 25 μ l of buffer). Place them back on ice.




2. Set up the Control:
 - Using a sterile transfer pipet, transfer 0.25 ml (250 μ l) of cell suspension from the tube "Cells" to the tube "Control Buffer".
 - Carefully place the pipet back into the wrapper.
 - Cap the tube; mix by tapping. Put the tube back on ice.



3. Set up the transformation:
 - Using the same pipet from Step 2, transfer 0.25 ml (250 μ l) of cell suspension from the tube "Cells" to the tube "pGAL DNA".
 - Cap the tube; mix by tapping. Put the tube back on ice.



4. Incubate the cells prepared in steps 1 - 3 on ice for **10 minutes**.



5. Place both transformation tubes at **42°C for 90 seconds**.

This heat shock step facilitates the entry of DNA in bacterial cells.

6. Return both tubes to the ice bucket and incubate for **1 minute**.




Transformation of *E. coli*

7. Add **0.75 ml** of the Recovery Broth to the tube "Control Buffer".

Add the recovery broth with a sterile 1 ml pipet. **Avoid touching the cells with the pipet.**

8. Add **0.75 ml** of the Recovery Broth to the tube "pGAL DNA".



Quick Reference:

DNA and competent cells are combined in a 0.25 ml suspension. After the cells have incubated with the DNA, growth medium (recovery broth) is added. Bacterial cells continue to grow through the recovery process, during which time the cell wall is repaired. Cells recover and begin to express the antibiotic resistance gene.


9. Incubate the closed tubes in a 37°C water bath for **30 minutes** for a recovery period.
10. While the tubes are incubating, label 3



agar plates as indicated below. Write on the bottom or side of the petri plate.

- Label one unstriped plate: X-GAL/Control 1
- Label one striped plate: AMP/X-GAL/Control 2
- Label one striped plate: AMP/X-GAL/pGAL
- Put your initials or group number on all the plates.

11. After the recovery period, remove the tubes from the water bath and place them on the lab bench. Proceed to plating the cells for



Duplication of this document, in conjunction with use of accompanying reagents, is permitted for classroom/laboratory use only. This document, or any part, may not be reproduced or distributed for any other purpose without the written consent of EDVOTEK, Inc. Copyright © 1989, 90, 93, 94, 97, 98, 2000, 2005, EDVOTEK, Inc., all rights reserved
EVT 005077K


Transformation of *E. coli*

PLATING THE CELLS

Plating cells from the tube labeled "Control":

12. Use a fresh, sterile 1 ml pipet to transfer recovered cells from the tube "Control Buffer" to the middle of the following plates:

- 0.25 ml to the plate labeled X-GAL/Control 1
- 0.25 ml to the plate labeled AMP/XGAL/Control 2



13. Spread the cells over the entire plate with a sterile inoculating loop.

14. Cover both control plates and allow the liquid to be absorbed.

To avoid contamination when plating, do not set the lid down on the lab bench -- Lift the lid of the plate only enough to allow spreading. Be careful to avoid gouging the loop into the agar.

Plating cells from the tube labeled "pGAL DNA"

15. Use a fresh, sterile 1 ml pipet to transfer recovered cells from the tube "pGAL DNA" to the middle of the following plate:

- 0.25 ml to the plate labeled AMP/X-GAL/pGAL

16. Spread the cells with a sterile inoculating loop.

17. Cover the plate and allow the liquid to be absorbed (approximately 15-20 minutes).

Reminder:
Follow proper procedures for disposal of contaminated materials.

Transformation of *E. coli***Reminder:**

Follow proper procedures for disposal of contaminated materials.

PREPARING PLATES FOR INCUBATION


20. Stack your group's set of plates on top of one another and tape them together. Put your initials or group number on the taped set of plates.

The plates should be left in the upright position to allow the cell suspension to be absorbed by the agar.

21. Place the set of plates in a safe place designated by your instructor.

22. After the cell suspension is absorbed by the agar (approximately 15 -20 minutes), you or your instructor will place the plates in the **inverted** position (agar side on top) in a 37°C bacterial incubation oven for overnight incubation (15-20 hours).

If the cells have not been absorbed into the medium, it is best to incubate the plates upright. The plates are inverted to prevent condensation on the lid, which could drip onto the culture and may interfere with experimental results.

**VIEWING PLATES AFTER INCUBATION**

23. Proceed to analyzing your results.
24. After analyzing your results, follow proper procedures for disposal of contaminated materials.



Duplication of this document, in conjunction with use of accompanying reagents, is permitted for classroom/laboratory use only. This document, or any part, may not be reproduced or distributed for any other purpose without the written consent of EDVOTEK, Inc. Copyright © 1989, 90, 93, 94, 97, 98, 2000, 2005, EDVOTEK, Inc., all rights reserved
EVT 005077K

Experiment Results and Analysis

LABORATORY NOTEBOOK RECORDINGS:

Address and record the following in your laboratory notebook or on a separate worksheet.

Before starting the Experiment:

- Write a hypothesis that reflects the experiment.
- Predict experimental outcomes.

During the Experiment:

- Record (draw) your observations, or photograph the results.

Following the Experiment:

- Formulate an explanation from the results.
- Determine what could be changed in the experiment if the experiment were repeated.
- Write a hypothesis that would reflect this change.

ANSWER THESE QUESTIONS BEFORE ANALYZING YOUR RESULTS.

1. On which plate(s) would you find only genetically transformed bacterial cells? Explain.
2. What is the purpose of the control plates? Explain the difference between each and why it is necessary to run each.
3. Why would one compare plates AMP/X-GAL and AMP/X-GAL/pGAL?

Continued

Experiment Results and Analysis

Data Collection

5. Observe the results you obtained on your transformation and control plates.

Transformation Plate: + DNA

- AMP/X-GAL/pGAL

Control Plates: - DNA

- X-GAL/ Control 1
- AMP/X-GAL/ Control 2

6. Draw and describe what you observe. For each of the plates, record the following:

- How much bacterial growth do you observe? Determine a count.
- What color are the bacteria?
- Why do different members of your class have different transformation efficiency values?
- If you did not get any results, what factors could be attributed to this fact?



Duplication of this document, in conjunction with use of accompanying reagents, is permitted for classroom/laboratory use only. This document, or any part, may not be reproduced or distributed for any other purpose without the written consent of EDVOTEK, Inc. Copyright © 1989, 90, 93, 94, 97, 98, 2000, 2005, EDVOTEK, Inc., all rights reserved
EVT 005077K

Experiment Results and Analysis

DETERMINATION OF TRANSFORMATION EFFICIENCY

Transformation efficiency is a quantitative determination of how many cells were transformed per 1 µg of plasmid DNA. In essence, it is an indicator of how well the transformation experiment worked.

You will calculate the transformation efficiency from the data you collect from your experiment.

1. Count the number of colonies on the plate with ampicillin that is labeled:

AMP/X-GAL/pGAL

A convenient method to keep track of counted colonies is to mark the colony with a lab marking pen on the outside of the plate.

2. Determine the transformation efficiency using the formula:

$$\frac{\text{Number of transformants}}{\mu\text{g of DNA}} \times \frac{\text{final vol at recovery (ml)}}{\text{vol plated (ml)}} = \frac{\text{Number of transformants}}{\text{per } \mu\text{g}}$$

Example:

Assume you observed 40 colonies:

$$\frac{40 \text{ transformants}}{0.025 \mu\text{g}} \times \frac{1.0 \text{ ml}}{0.25 \text{ ml}} = \frac{6400 \text{ (} 6.4 \times 10^3 \text{) transformants}}{\text{per } \mu\text{g}}$$

Quick Reference for Expt. 221

25 ng (0.025 µg) of DNA is used.

The final volume at recovery is 1.0 ml.

The volume plated is 0.25 ml.

Study Questions

Answer the following study questions in your laboratory notebook or on a separate worksheet.

1. Did you observe any satellite colonies? Why are the satellite, feeder colonies white?
2. Why did the competent cells which did not receive DNA (control) fail to grow on the plates containing ampicillin?
3. Why are there so many cells growing on the X-GAL plate? What color are they?
4. What evidence do you have that transformation was successful?
5. What are some reasons why transformation may be unsuccessful?



Duplication of this document, in conjunction with use of accompanying reagents, is permitted for classroom/laboratory use only. This document, or any part, may not be reproduced or distributed for any other purpose without the written consent of EDVOTEK, Inc. Copyright © 1989, 90, 93, 94, 97, 98, 2000, 2005, EDVOTEK, Inc., all rights reserved
EVT 005077K

Instructor's Guide

Notes to the Instructor:

IMPORTANT READ ME!

Transformation experiments contain antibiotics which are used for the selection of transformed bacteria. Students who have allergies to antibiotics such as penicillin, ampicillin, kanamycin or tetracycline should not participate in this experiment.

ORGANIZING AND IMPLEMENTING THE EXPERIMENT

Class size, length of laboratory sessions, and availability of equipment are factors which must be considered in the planning and the implementation of this experiment with your students.

The guidelines that are presented in this manual are based on ten laboratory groups consisting of two, or up to four students. The following are implementation guidelines, which can be adapted to fit your specific set of circumstances. If you do not find the answers to your questions in this section, a variety of resources are available at the EDVOTEK web site. In addition, Technical Service is available from 9:00 am to 6:00 pm, Eastern time zone. Call 1-800-EDVOTEK for help from our knowledgeable technical staff.

Online Ordering now available

www.edvotek.com

Visit our web site for information about EDVOTEK's complete line of experiments for biotechnology and biology education.

Technical Service Department

Mon - Fri
9:00 am to 6:00 pm ET

FAX: (301) 340-0582
web: www.edvotek.com
email: edvotek@aol.com

Please have the following information:

- The experiment number and title
- Kit Lot number on box or tube
- The literature version number (in lower right corner)
- Approximate purchase date

Day 1: (Prior to the Lab)

- Prepare agar plates
- Prepare LyphoCells™ (overnight incubation).
- Dispense the DNA and control buffer

Day 2: (Day of Lab Experiment)

- Equilibrate water baths at 37°C and 42°C; incubation oven at 37°C
- Students transform cells and plate for overnight incubation.

Day 3: (Day after Lab Experiment)

- Students observe transformants and controls
- Students calculate transformation efficiency
- Follow clean up and disposal procedures as outlined in the Laboratory Safety section.



EDVOTEK - The Biotechnology Education Company ®
1-800-EDVOTEK • www.edvotek.com
24-hour FAX: (301) 340-0582 • email: edvotek@aol.com



Notes to the Instructor:**NATIONAL CONTENT AND SKILL STANDARDS**

By performing this experiment, students will develop skills necessary to do scientific inquiry, learn new techniques using several types of biotechnology equipment, and will learn standard procedures used in transformation. Analysis of the experiments will provide students the means to transform an abstract concept into a concrete explanation.

APPROXIMATE TIME REQUIREMENTS

1. The experiment requires the reconstitution and incubation of LyphoCells at 34-37°C for 16-24 hours before the laboratory (overnight incubation). Plan accordingly. For optimal results, incubate LyphoCells at 34-37°C for 19 hours.
2. The agar plates can be prepared several days in advance and stored inverted (agar side on top) in the refrigerator. Preparation requires approximately 1 hour.
3. Dispensing the DNA and control buffer requires approximately 30 minutes. This can be done the day before the lab and stored in the refrigerator.
4. Competent cells must be dispensed just prior to the lab experiment. If tubes are already labeled, dispensing will require approximately 15 minutes.
5. Allow ample time for the equilibration of the water baths at 37°C and 42°C and a bacterial incubation oven at 37°C on the day of the experiment.
6. Each group will perform the transformation experiment and plate a set of three bacterial cells. These procedures require approximately 50 minutes.
7. Overnight incubation of plates is approximately 15-20 hours at 37°C. Colonies will also appear between 24 - 48 hours at room temperature.
8. Follow disposal procedures as outlined in the section regarding Laboratory Safety.

LABORATORY NOTEBOOKS

It is highly recommended that students maintain a laboratory notebook to formulate hypotheses and to record experimental procedures and results.

- EDVOTEK Cat. # 1401, Laboratory DataBook is recommended.
- Guidelines for keeping a laboratory notebook is available at the EDVOTEK web site.



Duplication of this document, in conjunction with use of accompanying reagents, is permitted for classroom/laboratory use only. This document, or any part, may not be reproduced or distributed for any other purpose without the written consent of EDVOTEK, Inc. Copyright © 1989, 90, 93, 94, 97, 98, 2000, 2005, EDVOTEK, Inc., all rights reserved
EVT 005077K

Pre-Lab Preparations



Wear Hot Gloves
and Goggles during all
steps involving heating.

POUR AGAR PLATES
(Prior to the Lab experiment)


- For optimal results, prepare plates two days prior to plating and set aside the plates inverted at room temperature.
- If they are poured more than two days before use, they should be stored inverted in the refrigerator. Remove the plates from the refrigerator and store inverted for two days at room temperature before use.

Heat the ReadyPour™ Medium

1. Equilibrate a water bath at 60°C for step 5 below.
2. Loosen, but do not remove, the cap on the ReadyPour medium bottle to allow for the venting of steam during heating.

Caution: Failure to loosen the cap prior to heating or microwaving may cause the ReadyPour medium bottle to break or explode.
3. Squeeze and vigorously shake the plastic bottle to break up the solid agar into chunks
4. Heat the bottle of ReadyPour medium by one of the methods outlined below. When completely melted, the amber-colored solution should appear free of small particles.
 - A. Microwave method:
 - Heat the bottle on High for two 30 second intervals.
 - Using a hot glove, swirl and heat on High for an additional 25 seconds, or until all the medium is dissolved.
 - Using a hot glove, occasionally swirl to expedite melting.
 - B. Hot plate or burner method:
 - Place the bottle in a beaker partially filled with water.
 - Heat the beaker to boiling over a hot plate or burner.
 - Using a hot glove, occasionally swirl to expedite melting.
5. Allow the melted ReadyPour medium to cool. Placing the bottle in a 60°C water bath will allow the agar to cool, while preventing it from prematurely solidifying.

When the ReadyPour™ medium reaches approximately 60°C, the bottle will be warm to the touch but not burning hot.



Pre-Lab Preparations

Add reagents to medium which has been cooled. Hot medium will cause reagents such as ampicillin to rapidly decompose.

Label ("Stripe") the Plates

6. Use a lab marker to "stripe" the sides of twenty (20) 60 x15 mm petri dishes. This will provide an easy method of differentiating between plates with ampicillin and plates without ampicillin.
 - Open one sleeve of 20 plates and stack the plates neatly.
 - Start the marker at the bottom of the stack and move the marker vertically to the top plate to "stripe" the sides of the 20 plates.
 - These plates will be used for medium with ampicillin.
 - Do not stripe the second sleeve of plates. These will be the control plates.

Pour the Plates (after the medium has cooled)

7. Thaw and add all of the X-Gal solution (Component E) to the molten and cooled ReadyPour medium. Recap bottle and swirl to mix.
8. Use a fresh 10 ml pipet and pipet pump to pour 10 unstriped plates, 5 ml each. (See Quick Reference: Pouring Agar Plates.) Save the pipet for step 11.

Quick Reference: Pouring Agar Plates

- Use a sterile 10 ml pipet with a pipet pump to transfer the designated volume of medium to each petri plate. Pipet carefully to avoid forming bubbles.
- Rock the petri plate back and forth to obtain full coverage.
- If the molten medium contains bubbles, they can be removed by passing a flame across the surface of the medium.
- Cover the petri plate and allow the medium to solidify.



Duplication of this document, in conjunction with use of accompanying reagents, is permitted for classroom/laboratory use only. This document, or any part, may not be reproduced or distributed for any other purpose without the written consent of EDVOTEK, Inc. Copyright © 1989, 90, 93, 94, 97, 98, 2000, 2005, EDVOTEK, Inc., all rights reserved
EVT 005077K

Pre-Lab Preparations

Reminder:

Follow proper procedures for disposal of contaminated materials.

9. Add the ampicillin powder (entire contents of tube D) to the remaining molten ReadyPour medium.
10. Recap the bottle; swirl to completely dissolve the ampicillin powder.
11. Use the pipet from step 8 to pour 20 striped plates, 5 ml each. Pour extra plates with any remaining medium.
12. Allow the agar to cool and resolidify.

Note: If plates will be used within two days, store in a sealable plastic bag so the plates will not dry out. Store at room temperature, inverted.

Summary of Poured Plates:

- 10 Agar plates with X-GAL
no stripe plates: 5 ml each - ReadyPour
medium without ampicillin
- 20 Agar plates with X-GAL and AMP
striped plates: 5 ml each - ReadyPour
medium with ampicillin

Pre-Lab Preparations

PREPARATION OF COMPETENT CELLS

Day before the Lab

1. Use a 10 ml sterile pipet to add 6 ml sterile cell reconstitution medium (Component F) to the vial of LyphoCells.
2. Replace the rubber stopper and cap. Mix by inverting until the freeze dried plug is dissolved.
3. Shake the cell suspension vigorously and incubate the vial at 34-37°C for 16 - 24 hours (overnight) in an incubation oven. For optimal results, incubate LyphoCells for 19 hours.

Day of the Lab

4. Completely thaw the competency induction solvent (G) and place on ice. (If there is a white precipitate in the bottle, warm it in a 37°C waterbath to dissolve the precipitate.)
5. Mix and resuspend the vial of incubated cells by inverting and gently shaking. Place the vial on ice for 10 minutes.
6. Use a 10 ml sterile pipet to add 3 ml of ice cold competency induction solvent (G) to the vial of cells.

The competency induction solvent is very viscous. Make sure that a portion of the solvent is not left on the walls of the pipet.

7. Mix the cells and induction solvent thoroughly by inverting the vial several times. The solution should have no dense layers, "streams" or globules (i.e. it should be a uniform suspension).
8. Keep the cells on ice for a minimum of 30 minutes.

Cells can be kept on ice for up to 3 hours.



Duplication of this document, in conjunction with use of accompanying reagents, is permitted for classroom/laboratory use only. This document, or any part, may not be reproduced or distributed for any other purpose without the written consent of EDVOTEK, Inc. Copyright © 1989, 90, 93, 94, 97, 98, 2000, 2005, EDVOTEK, Inc., all rights reserved
EVT 005077K

Pre-Lab Preparations

Keep cells on ice during dispensing procedures.

Dispensing the Cells just prior to the Experiment

11. Mix the cells by inversion to obtain an even suspension.
12. Use a sterile 1 ml pipet to aliquot 0.7 ml of cells to 10 ice cold tubes labeled "Cells"
13. Cap the tubes and keep them on ice.

PREPARATION OF DNA AND CONTROL BUFFER

1. Place the tubes of supercoiled pGAL DNA (Component B) and control buffer (Component C) on ice.
2. Before dispensing DNA and Control buffer, tap the tubes until all the sample is at the tapered bottom of the tube.

Summary of Reagent Preparations

		Dispense	
10 plates	Agar plates with X-GAL (No Stripe)	5 ml	
20 plates	Agar plates with AMP and X-GAL (Striped)	5 ml	
10 tubes	pGAL DNA	25 µl	on ice
10 tubes	Control Buffer	25 µl	on ice
10 tubes	Cells	0.7 ml	on ice

3. Using an automatic micropipet, dispense 25 µl of the supercoiled pGAL DNA to each of the 10 microtest tubes labeled "pGAL DNA".
4. Cap the tubes and place them on ice.
5. Using a FRESH micropipet tip, dispense 25 µl of control buffer to each of 10 microtest tubes labeled "Control Buffer".
6. Cap the tubes and place them on ice.

Each Group Requires:

1	tube	pGAL™ DNA	25 µl	on ice
1	tube	Cells	0.7 ml	on ice
1	tube	Control Buffer	25 µl	on ice
1	agar plate	- no stripe		
2	agar plates	- striped		
4	sterile 1 ml	pipets		
2	sterile inoculating	loops		
1	tube	Recovery broth (optional)	1.5 ml	
2	micropipet tips	(if using automatic micropipets)		

Recovery Broth (Optional)

Set up a classroom pipeting station, or dispense recovery broth (optional).

7. Dispense 1.5 ml Recovery Broth into 10 sterile tubes labeled "Recovery Broth" using a sterile pipet.
8. Cap the tubes and place them in the refrigerator if not to be used immediately.

Avoiding Common Pitfalls

1. When heating the ReadyPour medium, make sure it does not boil over and cause the volume to be reduced. Watch the bottle very carefully and remove it from heat if you see signs of the medium boiling over.
2. If the plates are made fresh, the plated cells will take longer to be absorbed into the medium. Invert the plates only after the cell suspension has been absorbed.
3. Do not discard the tubes containing transformed bacteria. After plating an aliquot on selection plates, set the tubes in a rack and leave on the lab bench overnight. If for some reason, transformants do not grow on the selection plates, the cell pellet can be plated as outlined below.



Technical Service Department

Mon - Fri
9:00 am to 6:00 pm ET

FAX: (301) 340-0582
web: www.edvotek.com
email: edvotek@aol.com

Please have the following information:

- The experiment number and title
- Kit Lot number on box or tube
- The literature version number (in lower right corner)
- Approximate purchase date

- Collect the bacterial cell pellet by centrifugation in a microcentrifuge. If a microcentrifuge is not available, let the bacteria collect by gravity and do not disturb.
- Remove all except 0.1 to 0.2 ml of medium (supernatant).
- Resuspend cell pellet in remaining medium.
- Spread entire contents of tube on selection medium.
- Incubate plate as before, 15 to 24 hours at 37°C.



Duplication of this document, in conjunction with use of accompanying reagents, is permitted for classroom/laboratory use only. This document, or any part, may not be reproduced or distributed for any other purpose without the written consent of EDVOTEK, Inc. Copyright © 1989, 90, 93, 94, 97, 98, 2000, 2005, EDVOTEK, Inc., all rights reserved
EVT 005077K

Experiment Results and Analysis

ANSWER THESE QUESTIONS BEFORE ANALYZING YOUR RESULTS.

1. On which plate(s) would you find only genetically transformed bacterial cells? Explain.

The bacteria growing on the plate labeled AMP/X-GAL/pGAL are transformed cells since only those cells that have taken up the plasmid which expresses the ampicillin resistance gene and the fluorescent gene(s) will survive on the plates which contain ampicillin.

2. What is the purpose of the control plates? Explain the difference between each and why it is necessary to run each.



Control plates help interpret the experimental results. There are two control plates in this experiment. The control plate that is labeled AMP/X-GAL/Control 2 shows that cells without the plasmid will not grow in the presence of ampicillin. The control plate X-GAL/Control 1 shows that host cells were not damaged during the transformation process and therefore are able to grow on agar plates that do not contain ampicillin.

3. Why would one compare plates AMP/X-GAL and AMP/X-GAL/pGAL?

Cells not treated with the plasmid will not grow on the plate with ampicillin (AMP/X-GAL) because they are not expressing the ampicillin resistance gene. However, cells treated with the plasmid will grow on the AMP/X-GAL/pGAL plate because they are expressing the ampicillin resistance gene.

Experiment Results and Analysis


IDEALIZED SCHEMATIC OF RESULTS.




Duplication of this document, in conjunction with use of accompanying reagents, is permitted for classroom/laboratory use only. This document, or any part, may not be reproduced or distributed for any other purpose without the written consent of EDVOTEK, Inc. Copyright © 1989, 90, 93, 94, 97, 98, 2000, 2005, EDVOTEK, Inc., all rights reserved
EVT 005077K

Study Questions and Answers


- 1. Did you observe any satellite colonies? Why are the satellite, feeder colonies white?**
Not all cells become transformed. Some of the untransformed cells can grow around the colonies containing transformed cells. The transformed cells secrete β -lactamase which clears the surrounding medium of ampicillin. Therefore, some of the untransformed cells grow. The satellites are white since they did not incorporate pGAL DNA which contains the gene that will allow the cell to have a functional β -galactosidase.
- 2. Why did the competent cells which did not receive DNA (control) fail to grow on the plates containing ampicillin?**
Without pGAL DNA, they are not ampicillin-resistant because they do not make β -lactamase.
- 3. Why are there so many cells growing on the X-GAL plate? What color are they?**
Cells have not been challenged with antibiotic and therefore virtually all survive. They are white in color because they do not have pGAL DNA, and therefore, do not have a functional β -galactosidase.
- 4. What evidence do you have that transformation was successful?**
A successful transformation will show colonies on the plate labeled AMP/X-GAL/pGAL. An unsuccessful transformation will not show any colonies.
- 5. What are some reasons why transformation may be unsuccessful?**
Unsuccessful transformations could be the result of many things, including: 1) not adding the plasmid to the host cells in the + pGAL DNA tube, or 2) not adding a colony of bacteria to the + pGAL DNA tube, and 3) improper timing of the heat shock step.


 <div> Material Safety Data Sheet May be used to comply with OSHA's Hazard Communication Standard. 29 CFR 1910.1200 Standard must be consulted for specific requirements. </div>			
IDENTITY (As Used on Label and List) Ampicillin		Note: Blank spaces are not permitted. If any item is not applicable, or no information is available, the space must be marked to indicate that.	
Section I			
Manufacturer's Name EDVOTEK, Inc. Address (Number, Street, City, State, Zip Code) 14676 Rothgeb Drive Rockville, MD 20850		Emergency Telephone Number (301) 251-5990 Telephone Number for information (301) 251-5990 Date Prepared 07/01/03 Signature of Preparer (optional)	
Section II - Hazardous Ingredients/Identify Information			
Hazardous Components [Specific Chemical Identity; Common Name(s)] Ampicillin		OSHA PEL ACGIH TLV Other Limits Recommended % (Optional) No data	
CAS# 7177-48-2			
Section III - Physical/Chemical Characteristics			
Boiling Point	No data	Specific Gravity (H ₂ O = 1)	No data
Vapor Pressure (mm Hg.)	No data	Melting Point	No data
Vapor Density (AIR = 1)	No data	Evaporation Rate (Butyl Acetate = 1)	No data
Solubility in Water Slightly soluble			
Appearance and Odor Odorless, white crystalline powder			
Section IV - Physical/Chemical Characteristics N.D. = No data			
Flash Point (Method Used) No data		Flammable Limits	LEL N.D. UEL N.D.
Extinguishing Media Dry chemical, carbon dioxide, water spray or regular foam			
Special Fire Fighting Procedures Move container from fire area if possible. Do not scatter spilled material with water streams.			
Unusual Fire and Explosion Hazards Avoid breathing vapors.			

Section V - Reactivity Data			
Stability	Unstable		Conditions to Avoid
	Stable	X	Incompatibles
Incompatibility Strong oxidizers			
Hazardous Decomposition or Byproducts Toxic oxides of carbon, nitrogen and sulfur			
Hazardous Polymerization	May Occur		Conditions to Avoid
	Will Not Occur	X	Incompatibles
Section VI - Health Hazard Data			
Route(s) of Entry: Inhalation? Yes Skin? Yes Ingestion? Yes			
Health Hazards (Acute and Chronic) Sensitizers may result in allergic reaction			
Carcinogenicity: No data NTP? IARC Monographs? OSHA Regulation?			
Signs and Symptoms of Exposure Repeated exposure may result in sensitization and possible anaphylactic shock.			
Medical Conditions Generally Aggravated by Exposure No data			
Emergency First Aid Procedures Ingestion: Allergic symptoms.			
Eyes/Skin: Flush with water Inhalation: Move to fresh air			
Section VII - Precautions for Safe Handling and Use			
Steps to be Taken in case Material is Released for Spilled Wear suitable protective clothing. Sweep up and place in suitable container for later disposal. Do not flush spilled material down sink.			
Waste Disposal Method Observe all federal, state, and local regulations			
Precautions to be Taken in Handling and Storing Keep away from incompatible substances			
Other Precautions None			
Section VIII - Control Measures			
Respiratory Protection (Specify Type)			
Ventilation	Local Exhaust	Yes	Special None
	Mechanical (General)	No	Other None
Protective Gloves	Yes	Eye Protection	Splash or dust proof
Other Protective Clothing or Equipment Eye wash			
Work/Hygienic Practices Wear protective clothing and equipment to prevent contact.			

 <div> Material Safety Data Sheet May be used to comply with OSHA's Hazard Communication Standard. 29 CFR 1910.1200 Standard must be consulted for specific requirements. </div>			
IDENTITY (As Used on Label and List) X-Gal in solvent		Note: Blank spaces are not permitted. If any item is not applicable, or no information is available, the space must be marked to indicate that.	
Section I			
Manufacturer's Name EDVOTEK, Inc. Address (Number, Street, City, State, Zip Code) 14676 Rothgeb Drive Rockville, MD 20850		Emergency Telephone Number (301) 251-5990 Telephone Number for information (301) 251-5990 Date Prepared 02/08/05 Signature of Preparer (optional)	
Section II - Hazardous Ingredients/Identify Information			
Hazardous Components [Specific Chemical Identity; Common Name(s)] This product contains no hazardous materials as defined by the OSHA Hazard Communication Standard.		OSHA PEL ACGIH TLV Other Limits Recommended % (Optional)	
Section III - Physical/Chemical Characteristics			
Boiling Point at 760mm Hg	189C	Specific Gravity (H ₂ O = 1)	No data
Vapor Pressure (mm Hg.) 20C	0.46	Melting Point	18C
Vapor Density (AIR = 1)	2.7	Evaporation Rate (Butyl Acetate = 1)	N/A
Solubility in Water Complete (100%)			
Appearance and Odor Colorless liquid, faint odor			
Section IV - Physical/Chemical Characteristics			
Flash Point (Method Used) (closed cup) 88C (192F)		Flammable Limits	LEL 3% 43% UEL
Extinguishing Media Water spray, carbon dioxide, dry chemical, ordinary foam			
Special Fire Fighting Procedures Wear SCBA with full facepiece operated in positive pressure mode. If possible, move container from fire area			
Unusual Fire and Explosion Hazards Vapors may flow along surfaces to distant ignition sources and flashback. Closed containers exposed to heat may explode. Contact w/ strong oxidizers may cause fire.			

Section V - Reactivity Data			
Stability	Unstable		Conditions to Avoid
	Stable	X	Heat, flame, other sources of ignition
Incompatibility Strong oxidizing agents, active halogen compounds, alkali metals			
Hazardous Decomposition or Byproducts Oxides of sulfur, mercaptane, active halogen compounds, alkali metals			
Hazardous Polymerization	May Occur		Conditions to Avoid
	Will Not Occur	X	None
Section VI - Health Hazard Data			
Route(s) of Entry: Inhalation? Yes Skin? Yes Ingestion? Yes			
Health Hazards (Acute and Chronic) Chronic effects: kidney damage, liver damage			
Carcinogenicity: None NTP? No IARC Monographs? No OSHA Regulation?			
Signs and Symptoms of Exposure Skin/eye irritation, headache, nausea, vomiting, dizziness			
Medical Conditions Generally Aggravated by Exposure Skin disorders			
Emergency First Aid Procedures Flush skin/eyes with large amounts of water. If ingested do not induce vomiting.			
Section VII - Precautions for Safe Handling and Use			
Steps to be Taken in case Material is Released for Spilled Shut off ignition sources. Wear protective clothing. Use water spray to reduce vapors. Take up with sand or other noncombustible absorbent material and dispose of properly			
Waste Disposal Method Dispose in accordance with all applicable federal, state, and local environmental regulations			
Precautions to be Taken in Handling and Storing Store in cool, dry, well ventilated flammable liquid storage area or cabinet, store above 20C			
Other Precautions Product may solidify at room temperature			
Section VIII - Control Measures			
Respiratory Protection (Specify Type) Chemical cartridge respirator w/ organic vapor cartridge.			
Ventilation	Local Exhaust	Yes	Special None
	Mechanical (General)	Yes	Other None
Protective Gloves	Butyl rubber gloves	Eye Protection	Splash Proof Safety goggles
Other Protective Clothing or Equipment Uniform or apron			
Work/Hygienic Practices Avoid contact with skin/eyes Keep container tightly closed.			

<div><div>Material Safety Data Sheet May be used to comply with OSHA's Hazard Communication Standard. 29 CFR 1910.1200 Standard must be consulted for specific requirements.</div></div>				Section V - Reactivity Data			
Stability		Unstable		Conditions to Avoid			
		Stable	X	Avoid incompatibles			
Incompatibility							
Strong oxidizers							
Hazardous Decomposition or Byproducts							
Sulfur dioxide, mercaptans, carbon monoxide, carbon dioxide, formaldehyde							
Hazardous Polymerization		May Occur		Conditions to Avoid			
		Will Not Occur	X				
Section VI - Health Hazard Data							
Route(s) of Entry:		Inhalation?		Skin?			
		Yes		Yes			
				Yes	Ingestion?		
Health Hazards (Acute and Chronic)		Inhalation/Ingestion: Nausea and vomiting					
		Skin/eye contact: Rapid absorption causing irritation					
Carcinogenicity:		NTP?		IARC Monographs?	OSHA Regulation?		
		None identified					
Signs and Symptoms of Exposure							
Irritation							
Medical Conditions Generally Aggravated by Exposure							
Skin disorders							
Emergency First Aid Procedures							
		Ingestion: Call medical help, do not induce vomiting		Inhalation: remove to fresh air			
		Skin/eye contact: Flush w/ water					
Section VII - Precautions for Safe Handling and Use							
Steps to be Taken in case Material is Released for Spilled							
Wear protective clothing. Take up with sand or other absorbant and place in container							
Dispose of properly.							
Waste Disposal Method							
Observe all federal, state, and local regulations.							
Precautions to be Taken in Handling and Storing							
Avoid contact							
Other Precautions							
None							
Section VIII - Control Measures							
Respiratory Protection (Specify Type)							
SCBA							
Ventilation		Local Exhaust	Yes	Special	None		
		Mechanical (General)	Yes	Other	None		
Protective Gloves		Butyl rubber gloves		Eye Protection			
				Safety goggles			
Other Protective Clothing or Equipment		Uniform or apron					
Work/Hygienic Practices		Avoid contact					

<div><div>Material Safety Data Sheet May be used to comply with OSHA's Hazard Communication Standard. 29 CFR 1910.1200 Standard must be consulted for specific requirements.</div></div>				
IDENTITY (As Used on Label and List)		Note: Blank spaces are not permitted. If any item is not applicable, or no information is available, the space must be marked to indicate that.		
Solvent for Induction of Competency				
Section I				
Manufacturer's Name		Emergency Telephone Number		
EDVOTEK, Inc.		(301) 251-5990		
Address (Number, Street, City, State, Zip Code)		Telephone Number for information		
14676 Rothgeb Drive		(301) 251-5990		
Rockville, MD 20850		Date Prepared		
		05-25-05		
		Signature of Preparer (optional)		
Section II - Hazardous Ingredients/Identify Information				
Hazardous Components [Specific Chemical Identity: Common Name(s)]		OSHA PEL	ACGIH TLV	Other Limits Recommended % (Optional)
None		-----	Not established	-----
Section III - Physical/Chemical Characteristics				
Boiling Point	No data	Specific Gravity (H ₂ O = 1)	No data	
Vapor Pressure (mm Hg.)	No data	Melting Point	No data	
Vapor Density (AIR = 1)	No data	Evaporation Rate (Butyl Acetate = 1)	No data	
Solubility in Water		Soluble		
Appearance and Odor		Clear liquid		
Section IV - Physical/Chemical Characteristics				
Flash Point (Method Used)	No data	Flammable Limits	LEL No data	UEL No data
Extinguishing Media				
Use water spray, alcohol foam, dry chemical, or carbon dioxide				
Special Fire Fighting Procedures				
Wear protective equipment and SCBA with full facepiece. Move container from fire area if possible.				
Unusual Fire and Explosion Hazards				
Vapors may flow along surfaces and flash back.				